

Mission®

Cholesterol Test Devices Package Insert

3-in-1 Lipid Panel	CHOL Total Cholesterol	TRIG Triglycerides	HDL High Density Lipoprotein
REF C131-2041	REF C131-2011	REF C131-2021	REF C131-2031
MODEL CCS-114	MODEL CCS-111	MODEL CCS-112	MODEL CCS-113

For testing cholesterol in human whole blood, plasma or serum.
For in vitro diagnostic use only.

INTENDED USE

The Mission® Cholesterol Test Devices work with the Mission® Cholesterol Meter to measure the lipid concentration in whole blood, plasma and serum during professional testing. For self-testing, use fingerip blood. The 3-in-1 Lipid Panel is used to measure the concentrations of Total Cholesterol (CHOL), High Density Lipoprotein (HDL) and Triglycerides (TRIG) simultaneously. It is also used to calculate LDL, CHOL/HDL and CHD risk. Note: CHD calculation function is only for professional use. Refer to the Mission® Cholesterol Monitoring System User's Manual for detailed instructions.

These separate test devices are also available, which can measure the concentrations of CHOL, HDL, and TRIG. Lipid measurements are used in the diagnosis and treatment of atherosclerotic coronary artery disease and in the diagnosis of metabolic disorders involving lipids and lipoproteins.

MEASUREMENT RANGE

Test Type	Measurement Range
Total Cholesterol	100-500 mg/dL (2.59-12.93 mmol/L)
High Density Lipoprotein	15-100 mg/dL (0.39-2.59 mmol/L)
Triglycerides	45-650 mg/dL (0.51-7.34 mmol/L)

For total cholesterol and high density lipoprotein, 1 mmol/L = 38.66 mg/dL; for triglycerides, 1 mmol/L = 88.6 mg/dL. Results below the ranges will show "<", and results above the ranges will show ">". When concentrations of specimens are above the test ranges, values for CHOL/HDL, LDL will display "...".

PRINCIPLE AND REFERENCE VALUES

Mission® Cholesterol Test Devices use a timed-endpoint method to measure the Total Cholesterol (CHOL), High Density Lipoprotein (HDL) and Triglyceride (TRIG) concentrations in whole blood, serum or plasma. The concentration of Low Density Lipoprotein (LDL) is calculated by the values of CHOL, TRIG and HDL. The system monitors the change in absorbance at 635 nm at a fixed time interval. The change in absorbance is directly proportional to the concentration of lipid in the specimen.

CHOL: In the reaction, cholesterol esterase hydrolyzes cholesterol esters to free cholesterol and fatty acids. The free cholesterol is oxidized to cholesterol-3-one and hydrogen peroxide by cholesterol oxidase. Peroxidase catalyzes the reaction of hydrogen peroxide with 4-aminoantipyrine and phenol to produce a colored quinoneimine product.

HDL: The dextran sulphate/Mg²⁺ on the test device precipitates the chylomicrons, VLDL and LDL, leaving HDL in the specimen. The cholesterol concentration of this HDL is then determined enzymatically, the same as CHOL.

TRIG: Triglycerides in the specimen are hydrolyzed to glycerol and free fatty acids by the action of lipase. A sequence of three coupled enzymatic steps using glycerol kinase (GK), glycerophosphate oxidase (GPO), and horseradish peroxidase (HPO) causes the oxidative coupling of 4-aminoantipyrine to form a blue dye.

LDL: When the concentration of TRIG in the specimen is equal to or lower than 400mg/dL, LDL concentration can be calculated by the meter with the following equation:
LDL = CHOL - HDL - TRIG/2 (mmol/L); LDL = CHOL - HDL - TRIG/5 (mg/dL)

Calculated LDL is an estimate of LDL.
Reference values are listed in the chart below^{2,4}:

Tests	Desirable	Borderline High	High
Total Cholesterol (CHOL)	<5.2 mmol/L (<200 mg/dL)	5.2-6.2 mmol/L (200-240 mg/dL)	>6.2 mmol/L (240 mg/dL)
High Density Lipoprotein (HDL)	Men: ≥1.5 mmol/L (60-40 mg/dL) Women: ≥1.3 mmol/L (60-50 mg/dL)	Men: 1.5-1.3 mmol/L (60-50 mg/dL) Women: <1.3 mmol/L (50 mg/dL)	Men: <1.0 mmol/L (40 mg/dL) Women: <1.3 mmol/L (50 mg/dL)
Triglycerides (TRIG)	<1.7 mmol/L (<150 mg/dL)	1.7-2.3 mmol/L (150-200 mg/dL)	>2.3 mmol/L (200 mg/dL)
Low Density Lipoprotein (LDL)	<3.4 mmol/L (<130 mg/dL)	3.4-4.1 mmol/L (130-160 mg/dL)	>4.1 mmol/L (160 mg/dL)

Reference ranges may vary between laboratories. Every laboratory should establish its own reference range as needed.¹ Blood lipid levels will have big physiological fluctuations depending on food consumed or exercise.

REAGENTS AND PERFORMANCE CHARACTERISTICS

Based on the dry weight at the time of impregnation, the concentrations given may vary within manufacturing tolerances.

Tests	Components
Total Cholesterol	Cholesterol esterase>0.3U; cholesterol oxidase>0.16U; POD(horseradish)>0.6U; ascorbate oxidase>0.6U;
High Density Lipoprotein	4-aminoantipyrine>0.06mg; Maos>0.06mg; buffer Magnesium chloride>0.1mg; dextran sulphate>0.6U; cholesterol oxidase>0.6U; 4-aminoantipyrine>0.06mg; Maos>0.06mg; buffer
Triglycerides	Liponitron lipase>0.35U; glycerol kinase>0.5U; glycerol phosphate oxidase>0.1U; POD(horseradish)>0.6U; ATP>0.2mg; ascorbate oxidase>0.5U 4-aminoantipyrine>0.09mg; Maos>0.06mg; buffer

The performance characteristics of these optical lipid devices have been determined in both laboratory and clinical tests. This test has been developed to be specific for the parameters to be measured with the exception of the interferences listed. Refer to the Limitations section for detailed information.

PRECAUTIONS

- For in vitro diagnostic use only.
- The test devices should remain in the original package until use.
- Do not use after the expiration date.
- Use the test device immediately after removing it from the foil pouch.
- Do not touch the reagent area of the test device.
- Discard any discolored or damaged test devices.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test device should be discarded according to local regulations after testing.
- Check the code chip before performing a test. Make sure to use the code chip that is included with the package of test devices. Insert the code chip into the code chip slot. The code chip slot is located on the left side of the meter.
- Check that the specimen type displayed on the meter LCD is the same as the specimen type tested. Whole blood samples have a two-digit test number that begins with the letter "B". Serum or plasma samples have the letter "S" in front of the test number. Decisions of medical relevance are not to be taken without consultation of a doctor. Changes to treatment should only be made after proper training.

Total precision is listed below:

Total Precision	Level I (n=60)	Level II (n=60)	Level III (n=60)
Mean (mg/dL)	90	192	321
SD (mg/dL) or %CV	3.89	3.20%	3.60%

Accuracy

The Cholesterol Test Devices were used by a trained technician to test heparin-preserved venous whole blood specimens from 78 participants. The same specimens were analyzed using a reference method (X). The results are compared below:

Specimen	Slope	Intercept	R	N
Venous whole blood	1.0243	-2.7846	0.994	78

High Density Lipoprotein

Specimen	Slope	Intercept	R	N
Venous whole blood	0.9728	1.6124	0.991	78

Triglycerides

Specimen	Slope	Intercept	R	N
Venous whole blood	0.9891	1.4849	0.993	78

In another study, heparinized venous whole blood, serum and heparinized plasma were collected from each patient and tested using a Cholesterol Test Device by a trained technician. A total of 40 patients took part in this study and results compared to those tested on the serum from same patients by the Abel-Kendall method (For CHOL) and DCM method (For HDL) in a Cholesterol Reference Method Laboratory Network (CRMLN) laboratory. The results are listed below.

Total Cholesterol

Specimen	Slope	Intercept	R	N
Venous whole blood	1.0286	-6.5223	0.998	40
Plasma	1.0336	-4.4486	0.998	40
Serum	1.0402	-6.145	0.999	40

High Density Lipoprotein

Specimen	Slope	Intercept	R	N
Venous whole blood	1.0334	-0.6386	0.995	40
Plasma	1.0441	-0.7255	0.995	40
Serum	1.0438	-0.8096	0.995	40

QUALITY CONTROL

For best results, performance of test devices should be confirmed by testing known specimens/controls whenever a new test is performed or whenever a new package is first opened. Each laboratory should establish its own goals for adequate standards of performance. Contact your local distributor for information on specific controls for this product. The following substances do not interfere with test results:

Substance	Amount	Substance	Amount
Acetaminophen	1324 µmol/L (20 mg/dL)	Cholesterol	12.9 mmol/L (500 mg/dL)
Ascorbic Acid	568 µmol/L (10 mg/dL)	Triglyceride	7.3 mmol/L (650 mg/dL)
Conjugated Bilirubin	240 µmol/L (20 mg/dL)	Uric Acid	0.6 mmol/L (10 mg/dL)
Creatinine	442 µmol/L (5 mg/dL)	Hemoglobin	2 g/L (200 mg/dL)
Ibuprofen	2425 µmol/L (50 mg/dL)	Dopamine	5.87 µmol/L (0.09 mg/dL)
Methylopa	71 µmol/L (1.5 mg/dL)	N-Acetylbenzoinoneimine	1.34 µmol/L (20mg/dL)

High concentrations of uric acid and ascorbic acid can lead to low measurements. Anticoagulants, such as heparin and EDTA, are recommended for use with venous whole blood. Do not use EDTA plasma, which leads to higher results. Do not use other anticoagulants, such as sodium citrate, sodium oxalate, sodium heparin, sodium fluoride, or sodium fluoride. Anticoagulant therapy may lower results. Venous blood is not recommended for use. Hemolyzed blood or blood from patients on thrombolytic therapy may lower results. Venous occlusion may increase results, and it is not recommended that blood be drawn from sites where veins are occluded.

BIBLIOGRAPHY

- Henry, J. B. Clinical Diagnosis and Management by Laboratory Methods. 15-290, 2001.
- Friedewald et al. Clin Chem. 1972. 18(6): 499-502
- National Cholesterol Education Program 2001 Guidelines, National Institutes of Health, May 2001.
- ATP III NCEP Guidelines for CHD Risk. JAMA. 2001. 285:2486-2509

INDEX OF SYMBOLS

Symbol	Consult instructions for use	Use-by date	Temperature limit
IVD	In vitro diagnostic medical device	Batch code	Control range
CODE	Code number	Manufacturer	Catalogue number
REF	Contains sufficient for <N> tests	Model number	Do not reuse
EC REP	Authorized representative in the European Community		

ACOM

ACOM Laboratories, Inc.
5650 Oberlin Drive, #340
San Diego, CA 92121, USA

CE 0123

MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany

Number: 1151078001
Effective date: 2020-10-16